



DELUXE EXOFLEX ORDER FORM

175 Park Street
Livermore Falls, ME 04254

PH: (207) 897-5558
FX: (207) 897-1117

www.pinetreeorthopedic.com

Doctor/Facility: _____ Date: ___/___/___

Address: _____ Customer Ref # _____

Phone: (____) _____ - _____

Patient Name: _____

Left Right Bi - Lateral

Diagnosis/
Observations _____

Age: _____ Weight: _____ Male Female Shoe Size: _____ Shoe Style: _____

Brace Options

Please complete this entire section. Missing information may cause production delays.

Control: Medium (Semi Rigid) Maximum (Firm)
** Patient weight is used to determine material thickness. If no weight is given, material thickness will not be adjusted.

Height: 7" 9" As Marked on Cast
** Height is measured from base of heel to top of brace. **Cast must be 1" above desired brace height.**

Foot Plate Trim Length: Proximal to Met Heads Sulcus Full Foot Extension As Marked
**Length of Plastic Only

Padding Extension: Sulcus Full Foot Extension Cut to shoe size
**Leave blank for met length

Padding Thickness 1/8" Swirl EVA 1/16" Swirl EVA 1/8" Plastazote

Special Options Use Extrinsic Rearfoot Post (Additional charges apply)

Special Modifications/ Instructions

Cast Preparation

Ankle Correction: As Casted Correct to 90°
Forefoot Correction: As Casted Correct to Neutral

STS Sock Return

Yes, return a Medium Large STS Sock with this order.